

# Foot Restraints

The neoprene ankle straps are available in two versions. This standard version with velcro fasteners is used where a simple correction is required. It should also be used where uncontrolled actions could lead to injury. The strap would then open automatically when overloaded. For example when used by clients who suffer from extension spasms.

## Neoprene Ankle Straps



Art. No.	E 4 1 400=1	E 4 1 400=2	E 4 1 400=3
	Size 1	Size 2	Size 3
smallest ankle circumference	20 cm	26 cm	28 cm
biggest ankle circumference	26 cm	32 cm	34 cm
strap width	4,5 cm	5,5 cm	6,5 cm



## Footshells, plastic

These light yet sturdy foot shells are vacuum formed. The profile ensures stability while incorporating an attractive design. They may be used in conjunction with velcro, buckles and straps or toothed belts with ratchet fastening.



Description	Art. No.	Size*
Footshells (plastic) with straps, size 1	E 4 1 100=1	15 x 6 x 8 cm
Footshells (plastic) with straps, size 2	E 4 1 100=2	19 x 7 x 9 cm
Footshells (plastic) with straps, size 3	E 4 1 100=3	22 x 8 x 10 cm
Footshells (plastic) with straps, size 4	E 4 1 100=4	25 x 10 x 12 cm
Footshells (plastic) without straps, size 1	E 4 1 200=1	15 x 6 x 8 cm
Footshells (plastic) without straps, size 2	E 4 1 200=2	19 x 7 x 9 cm
Footshells (plastic) without straps, size 3	E 4 1 200=3	22 x 8 x 10 cm
Footshells (plastic) without straps, size 4	E 4 1 200=4	25 x 10 x 12 cm

## Footshells, aluminium



Description	Art. No.	Size*
Footshells (aluminium) with straps, size 1	E 4 1 210=1	15 x 6 x 8 cm
Footshells (aluminium) with straps, size 2	E 4 1 210=2	19 x 7 x 9 cm
Footshells (aluminium) with straps, size 3	E 4 1 210=3	22 x 8 x 10 cm
Footshells (aluminium) without straps, size 1	E 4 1 220=1	15 x 6 x 8 cm
Footshells (aluminium) without straps, size 2	E 4 1 220=2	19 x 7 x 9 cm
Footshells (aluminium) without straps, size 3	E 4 1 220=3	22 x 8 x 10 cm

The aluminium foot shells are extremely robust. They are powder coated and easy to clean. They can also be easily adapted for customised applications.

\*These given measurements are inside measurements: inside length x heel width x forefoot width

# Foot Restraints

## Flexi Footshells



<b>Description</b>	<b>Art. No.</b>	<b>Size*</b>
Flexi Footshells, size 1	E 4 1 250=1	20 x 7 x 9 cm
Flexi Footshells, size 2	E 4 1 250=2	26 x 9 x 11 cm
Flexi Footshells, size 3	E 4 1 250=3	30 x 12 x 13 cm

\*These given measurements are inside measurements: inside length x heel width x forefoot width