



Lap Belts

Lap belts should be used for correct pelvic positioning. They are not safety belts and should only be used for transportation in combination with approved restraint systems. These lap belts provide a wide range of possibilities for pelvic positioning. For children we suggest you use 25 – 40 mm lap belts and for adults 40 – 50 mm lap belts. For clients who display a high tonus the lap belt with metal buckle should be considered. We offer a larger range of standard and also custom solutions for almost all cases. Please contact us if you require a version made to your specifications.




Lap Belt with metal buckle

Art. No.	A	B
E 1 5 100=40	650 mm	40 mm




Lap Belt with plastic buckle

Art. No.	A	B
E 1 5 400=40	800 mm	40 mm
E 1 5 400=50	800 mm	50 mm




Neoprene Lap Belt with plastic buckle




Art. No.	A	B	C	D
E 1 5 410=25	100 mm	50 mm	25 mm	750 mm
E 1 5 410=30	125 mm	60 mm	30 mm	750 mm
E 1 5 410=40	150 mm	70 mm	40 mm	750 mm
E 1 5 410=50	170 mm	85 mm	50 mm	750 mm

Neoprene Lap Belt with plastic buckle and bandalier




Art. No.	A	B	C	D
E 1 5 411=25	100 mm	50 mm	25 mm	750 mm
E 1 5 411=30	125 mm	60 mm	30 mm	750 mm
E 1 5 411=40	150 mm	70 mm	40 mm	750 mm
E 1 5 411=50	170 mm	85 mm	50 mm	750 mm

Neoprene Lap Belt with metal buckle




Art. No.	A	B	C	D
E 1 5 412=25	100 mm	50 mm	25 mm	750 mm
E 1 5 412=38	150 mm	70 mm	40 mm	750 mm
E 1 5 412=48	170 mm	85 mm	50 mm	750 mm

Lap Belts



Neoprene Lap Belt with two plug-type buckles



Art. No.	A	B	C	D
E 1 5 420=25	180 mm	60 mm	25 mm	700 mm
E 1 5 420=30	260 mm	70 mm	30 mm	700 mm
E 1 5 420=40	300 mm	85 mm	40 mm	700 mm
E 1 5 420=50	350 mm	85 mm	50 mm	700 mm

Neoprene Lap Belt with two plug-type buckles and bandalier




Art. No.	A	B	C	D
E 1 5 421=25	180 mm	60 mm	25 mm	700 mm
E 1 5 421=30	260 mm	70 mm	30 mm	700 mm
E 1 5 421=40	300 mm	85 mm	40 mm	700 mm
E 1 5 421=50	350 mm	85 mm	50 mm	700 mm


Neoprene Lap Belt with two plug-type buckles, bandalier and one side pull



Art. No.	A	B	C	D
E 1 5 422=30	260 mm	70 mm	30 mm	700 mm
E 1 5 422=40	300 mm	85 mm	40 mm	700 mm
E 1 5 422=50	350 mm	85 mm	50 mm	700 mm




4-Point Neoprene Lap Belt with plastic buckle



Art. No.	A	B	C	D	E	F
E 1 5 440=25	100 mm	50 mm	25 mm	750 mm	20 mm	650 mm
E 1 5 440=30	125 mm	60 mm	30 mm	750 mm	25 mm	650 mm
E 1 5 440=40	150 mm	70 mm	40 mm	750 mm	25 mm	650 mm
E 1 5 440=50	170 mm	85 mm	50 mm	750 mm	25 mm	650 mm

4-Point Neoprene Lap Belt with plastic buckle and bandalier



Art. No.	A	B	C	D	E	F
E 1 5 441=25	100 mm	50 mm	25 mm	750 mm	20 mm	650 mm
E 1 5 441=30	125 mm	60 mm	30 mm	750 mm	25 mm	650 mm
E 1 5 441=40	150 mm	70 mm	40 mm	750 mm	25 mm	650 mm
E 1 5 441=50	170 mm	85 mm	50 mm	750 mm	25 mm	650 mm