

Vacuum Bags

Our pre-formed and flat vacuum bags are available in four sizes (vacuum mattresses in two sizes). They are manufactured using a **special, latex-free material**. The anatomical form of the vacuum bags helps to simplify the casting process and to create a precise, individual cast of the client according to the therapeutic/ medical requirements.

The vacuum bags are filled with a fine dimensionally stable polystyrene granulate. A vacuum pump (we supply both manual and electric pumps from our stock) is used to evacuate the air out of the bags to create a vacuum, this results in a *firm cast of the client*. This process can be repeated by simply adding air to the bag, repositioning the client, until an optimal position has been achieved. Once the ultimate cast has been achieved the bag may now be used for a trial seating test in order to assess the clients positioning. Although the cast is rigid the anatomic form ensures that the client is comfortable during this short test phase. These vacuum bags are not suitable for long-term positioning or for permanent use. They are a tool used to assess a client and to make a cast only.

All vacuum bags are equipped with quick fitting couplings. Simply push both couplings together until an audible „click“ is heard. You have now made a sealed connection. To separate both couplings, simply press the release button and the two couplings will part.

Our Jumbo Holdall bags are available in two sizes and offer the best protection against damage. They are suitable for storage and transporting the bags.



Article No.	Description	Width
E 1 6 100=30R	Vacuum bag (latex free), light green back	30 cm
E 1 6 100=30S	Vacuum bag (latex free), light green seat	30 cm
E 1 6 100=35R	Vacuum bag (latex free), light green back	35 cm
E 1 6 100=35S	Vacuum bag (latex free), light green seat	35 cm
E 1 6 100=40R	Vacuum bag (latex free), light green back	40 cm
E 1 6 100=40S	Vacuum bag (latex free), light green seat	40 cm
E 1 6 100=45R	Vacuum bag (latex free), light green back	45 cm
E 1 6 100=45S	Vacuum bag (latex free), light green seat	45 cm



Article No.	Description	Measurements
E 1 6 140	Vacuum bag (latex free), light green flat	40 x 40 cm
E 1 6 145	Vacuum bag (latex free), light green flat	45 x 45 cm
E 1 6 170	Vacuum bag (latex free), light green flat	70 x 30 cm
E 1 6 200	Vacuum bag (latex free), light green flat	90 x 80 cm
E 1 6 220	Vacuum bag (latex free), light green flat	38 x 19 x 5 cm
E 1 6 300	Vacuum bag mattress (latex free)	190 x 90 cm
E 1 6 400	Vacuum bag mattress (latex free)	160 x 90 cm

