

# Vacuum Bags

Our pre-formed and flat vacuum bags are available in four sizes (vacuum mattresses in two sizes). They are manufactured using a **special, latex-free material**. The anatomical form of the vacuum bags helps to simplify the casting process and to create a precise, individual cast of the client according to the therapeutic/ medical requirements.

The vacuum bags are filled with a fine dimensionally stable polystyrene granulate. A vacuum pump (we supply both manual and electric pumps from our stock) is used to evacuate the air out of the bags to create a vacuum, this results in a *firm cast of the client*. This process can be repeated by simply adding air to the bag, repositioning the client, until an optimal position has been achieved. Once the ultimate cast has been achieved the bag may now be used for a trial seating test in order to assess the clients positioning. Although the cast is rigid the anatomic form ensures that the client is comfortable during this short test phase. These vacuum bags are not suitable for long-term positioning or for permanent use.

All vacuum bags are equipped with quick fitting couplings. Simply push both couplings together until an audible „click“ is heard. You have now made a sealed connection. To separate both couplings, simply press the release button and the two couplings will part.

Our Jumbo Holdall bags are available in two sizes and offer the best protection against damage. They are suitable for storage and transporting the bags.



Article No.	Description	Width
E 1 6 102=30R	Vacuum bag (latex free), blue back	30 cm
E 1 6 102=30S	Vacuum bag (latex free), blue seat	30 cm
E 1 6 102=35R	Vacuum bag (latex free), blue back	35 cm
E 1 6 102=35S	Vacuum bag (latex free), blue seat	35 cm
E 1 6 102=40R	Vacuum bag (latex free), blue back	40 cm
E 1 6 102=40S	Vacuum bag (latex free), blue seat	40 cm
E 1 6 102=45R	Vacuum bag (latex free), blue back	45 cm
E 1 6 102=45S	Vacuum bag (latex free), blue seat	45 cm



Article No.	Description	Measurements
E 1 6 140	Vacuum bag (latex free), blue flat	40 x 40 cm
E 1 6 145	Vacuum bag (latex free), blue flat	45 x 45 cm
E 1 6 202	Vacuum bag (latex free), blue flat	90 x 80 cm
E 1 6 222	Vacuum bag (latex free), blue flat	38 x 19 x 5 cm
E 1 6 302	Vacuum bag mattress (latex free)	190 x 90 cm
E 1 6 402	Vacuum bag mattress (latex free)	160 x 90 cm



# Vacuum Bags

Manual and electric pumps are offered as accessories. Both are fitted with tubes and quick coupling valves, so that you can modify your existing vacuum bags. Therefore, it is possible to adapt already existing equipment from other manufacturers to a uniform system.



Art. No.	Description	Measurements
E 1 6 020	Pump, manual with valve	
E 1 6 040	Pump, electric with valve	
E 1 6 099	Jumbo Holdall for vacuum bags	70 x 60 x 18 cm
E 1 6 399	Jumbo Holdall for vacuum mattresses	100 x 100 x 35 cm



Art. No.	Description
E 1 6 010	Valve (for bag)
E 1 6 011	Valve (for pump)

